CE 0483

EN

INSTRUCTIONS FOR USE

Description

Synotec[®] molecular matrix (sodium hyaluronate) is a 2.5% viscoelastic, sterile and clear single dose solution for injection into the joint cavity. The highly purified sodium hyaluronate is dissolved in a physiological buffer at pH 6.8 - 7.6.

Composition per syringe

Each syringe of 3 ml (Synotec[®] molecular matrix) contains:

75 mg

q.s.

q.s.

q.s.

Intended use

Treatment of the moderate osteoarthritis in knee, hip, shoulder, ankle and temporomandibular joint. Synotec[®] molecular matrix can also give extra viscoelastic support to other joints with synovial capsule.

Characteristics and mode of action

Hyaluronic acid is a native substance and hence very biocompatible. The most important properties of hyaluronic acid are protecting, lubricating and supporting of delicate cells and tissues.

Synotec[®] molecular matrix, containing cross-linked sodium hyaluronate, free from avian protein, may ameliorate the viscoelastic properties of synovial fluid, thus the lubrication and shock-absorbing properties of synovial fluid improved and the mechanical overload of the joints is reduced.

Dose and administration

Synotec[®] molecular matrix is injected into the affected joint only once. Several joints may be treated at the same time. Repeated treatment cycles may be administered as required at intervals of not less than 6 months. Take the pre-filled syringe out of the blister. Before use allow the preparation to equilibrate to room temperature for 20 minutes. Remove the tip cap. Mount the cannula by using the thumb and forefinger to hold firmly around both the glass syringe barrel and the Luer-lock adapter. Grasp the cannula shield with the other hand. To facilitate proper assembly, both push and rotate firmly. Suitable cannula size (e.g., 18-22 G). Work especially to open the cannula and to remove the cap on the syringe.

Contraindications and precautions

- Do not use in patients with ancertated hypersensitivity to the constituents in Synotec[®] molecular matrix.

- Synotec[®] molecular matrix must not be used when patients have septic arthritis. It is not indicated for use in other types of arthritis not related to osteoarthritis.

- Synotec[®] molecular matrix must only be injected by qualified medical personnel (orthopaedist, rheumatologist or equivalent). The dosage shall be prescribed by the physician, who will adjust the dosage specifically for the individual patient.

- The syringe is for single use only. Do not reuse.

- Reusable synovial fluid before injecting in the joint. Synotec[®] molecular matrix must be injected aseptically into the joint cavity. Avoid injections into blood vessels and surrounding tissues.

- General precautions for intra-articular injections shall be observed.

- Limited experience is available for the use of hyaluronic acid in children, pregnant and lactating women. Consult your doctor concerning Synotec[®] molecular matrix in these cases.

- Synotec[®] molecular matrix may contain traces of BDEs, which in rare cases can cause tissue irritation.

Side effects

Local secondary phenomena, such as pain, feeling of heat, redness of the skin (flush) and swelling may occur in a joint treated with Synotec[®] molecular matrix.

Interactions with other products

No inject anesthetics or other pharmaceutical preparations in the joint during the period of time the patient is treated with Synotec[®] molecular matrix.

Please note

The solution for injection in the syringe is sterile. For storage and transportation the syringe is protected by a blister package in a carton. The outer surface of the pre-filled syringe is clean, but not sterile.

- Do not use syringe, if the pre-filled syringe or blister cover is damaged. Strict aseptic administration technique must be followed.

Store between +2°C and +25°C. Protect from freezing and direct light.

Shelf life: see labeling/marketing on syringe and product carton.

Return to a physician only.

Store away from children!

Package size

1 Pre-filled syringe, containing 3 ml (Synotec[®] molecular matrix) of sterile solution for injection.

Do not re-use Steam sterilized Do not use if package is damaged Keep away from sunlight Use-by date

Consult instruction for use Store between +2°C and +25°C LOT Batch code

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